

SEE THE DIFFERENCE WITH THE

# Baylis RF *RadiOpaque*<sup>™</sup> Cannula

The *RadiOpaque*<sup>™</sup> cannula provides better visualization during positioning, resulting in a better targeted treatment.



Platinum adds value

# RadiOpaque™ Cannula

The *RadiOpaque™* cannula features a platinum band located at the distal end of the active tip. This defines the actual lesion location under fluoroscopy. The *RadiOpaque™* cannula provides better visualization during positioning, resulting in a better targeted treatment.

22 GAUGE ■  
21 GAUGE ■  
20 GAUGE ■  
18 GAUGE ■



## RadiOpaque™ Cannula


- 1 Enhanced visualization during fluoroscopy
- 2 Accurate positioning
- 3 Available in several size/tip configurations
- 4 Color-coded for quick reference
- 5 Disposable to prevent cross-contamination
- 6 Integral part of Baylis Medical RF Pain Management System.

Stylet

Platinum Band

Active Tip

Insulation



**Baylis**  
MEDICAL

Baylis Medical Company Inc.  
5959 Trans-Canada Highway  
Montreal, QC, Canada, H4T 1A1  
Tel. : (514) 488-9801  
Fax : (514) 488-7209  
info@baylismedical.com  
www.baylismedical.com

PM1004 Rev 08/08

Baylis Medical is a world leader in the development and distribution of high-technology pain management products. We conceive, innovate, develop, manufacture, and sell state-of-the-art medical products used around the world. Baylis Medical operates under the ISO 13485:2003 approved quality system.

© Copyright Baylis Medical Company Inc., 2008. The BMC logo and *RadiOpaque™* are trademarks or registered trademarks of Baylis Medical Company Inc. in the United States and/or other countries. Baylis Medical Company Inc. reserves the right to change specifications or to incorporate design changes without notice and without incurring any obligation relating to equipment previously manufactured or delivered.

Caution : Federal Law (USA) restricts the sale of these devices to or by the order of a physician.

Patents Pending and/or issued